


**1. Part Number:** M1010B-150

**2. Photograph:**



### 3. Product Information:

Description:	Metal Detectable V-Spatula 150mm (Type1) - Blue
Material of Construction:	Polypropylene (PP) <ul style="list-style-type: none"> <li>• Metal detectable grade</li> <li>• Only virgin material is used</li> <li>• PP conforms to FDA 21 CFR 177.1520</li> <li>• Blue colourant conforms to FDA 21 CFR 178.3297</li> <li>• Filler conforms to FDA 21 CFR 186.1374</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Injection moulded as a single piece
Overall Length:	248mm
Blade Length:	150mm
Blade Width:	20mm
Nominal Weight of Sampler:	19g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Number of samplers per Pack:	4
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	Not Applicable – reusable product
ATEX:	This device has not been ATEX tested

**4. Country of Origin:**

Made and packed in the UK.

**5. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Aug-24