

1. Part Number: M1016B-200

2. Photograph:



3. Product Information:

| | |
|---------------------------------|---|
| Description: | Metal Detectable Pallet Knife 200mm - Blue |
| Material of Construction: | Polypropylene (PP) <ul style="list-style-type: none"> • Metal detectable grade • Only virgin material is used • PP conforms to FDA 21 CFR 177.1520 • Blue colourant conforms to FDA 21 CFR 178.3297 • Filler conforms to FDA 21 CFR 186.1374 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free |
| Method of Construction: | Injection moulded as a single piece |
| Overall Length: | 310mm |
| Blade Length: | 200mm |
| Blade Width: | 45mm |
| Nominal Weight: | 92g |
| Moulding Environment: | Class 100,000 Cleanroom |
| Packing Environment: | Class 10,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Number of scoops per Pack: | 4 |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | Not Applicable – reusable product |
| ATEX: | This device has not been ATEX tested |

4. Country of Origin:

Made and packed in the UK.

5. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Aug-24