

1. Part Number: M1036B-20

2. Photograph:



3. Product Information:

Description:	Metal Detectable Long Handled Spoon 20ml - Blue
Material of Construction:	Polypropylene (PP) • Metal detectable grade • Only virgin material is used • PP conforms to FDA 21 CFR 177.1520 • Blue colourant conforms to FDA 21 CFR 178.3297 • Filler conforms to FDA 21 CFR 186.1374 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Injection moulded as a single piece
Nominal Bowl Volume:	20ml
Handle Length:	285mm
Bowl Length:	88mm
Bowl Width:	38mm
Nominal Weight:	34g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Number of scoops per Pack:	4
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	Not Applicable – reusable product
ATEX:	This device has not been ATEX tested

4. Country of Origin:

Made and packed in the UK.

5. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 22-Aug-24